

Bone Graft Substitutes

OrthoNovis BTB Morsels



Rapid Mineralization, Complete Absorption.

OrthoNovis BTB Morsels are bone graft substitutes designed with a continuous porous structure. The open structure and interconnected porosity facilitate rapid mineralization, leading to complete absorption and replacement by new bone during the healing process.

Properties:

- Continuous open porous structure enables rapid mineralization
- Available in a sealed cup or pre-loaded syringe with a Luer Lock attachment, facilitating the easy expression of blood directly to the morsels for improved handling during implantation
- Provides a bone void filler that resorbs and is replaced with bone during the healing process

Indications:

- Indicated only for bony voids or gaps that are not intrinsic to the stability of the bony structure
- Indicated to be gently packed into bony voids or gaps of the skeletal system. These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone

OrthoNovis BTB Morsels

5.0cc Morsels 1-2mm Cup NB1805NV

10cc Morsels 1-2mm Cup NB1810NV